

AutoScript TEM 1.18.0

Release Notes

PN 1555424

Revision A • June 2026 • Original Instructions

Confidential, Limited Rights

Copyright and Trademarks

Technical Publications

Technical Publications Team - Hillsboro

Copyright © 2026 by FEI Company, a part of Thermo Fisher Scientific. The information and materials contained herein are proprietary to Thermo Fisher and are provided for your organization's internal use on a need-to-know basis. They cannot be duplicated, published or disseminated to any third party without the express written consent of Thermo Fisher.

Trademark Acknowledgments

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. FEI and the FEI logo are registered trademarks of FEI Company (part of Thermo Fisher Scientific and its affiliates). All other trademarks belong to their respective owners.

Excel and Microsoft are registered trademarks of Microsoft Corporation.

Limited Rights

The following notice applies to the U.S. Government and other purchases with federal funds:

Contractor Name: Thermo Fisher Scientific

Contractor Address: 5350 NE Dawson Creek Drive, Hillsboro OR 97124

The Government's rights to use, modify, reproduce, release, perform, display, or disclose these technical data are restricted to those rights specified in DFARS 252.227-7015(b)(2), FAR 52.227-14(g)(2)(Alternate II) and FAR 12.211. Any reproduction of technical data or portions thereof marked with this legend must also reproduce the markings. Any person, other than the Government, who has been provided access to such data, must promptly notify the above named Contractor.

Document History

Rev A, June 2026

Introduction

Purpose

This document describes the Thermo Fisher Scientific AutoScript TEM software releases.

Audience and Scope

These release notes are intended for users of the Thermo Fisher Scientific AutoScript TEM software and those who manage the installation of Thermo Fisher Scientific AutoScript TEM software on the microscope computer.

This document describes the content of the most recent AutoScript TEM releases.

Hardware Requirements

The AutoScript TEM Server software must be installed on the Microscope PC.

The AutoScript TEM Client software can be installed on either the Microscope PC, or on the Support PC or Network PC

For successful installation and use, the AutoScript TEM Client software requires the following hardware specifications:

Specification	Minimum	Recommended
CPU	Intel Core i5 with at least 4 cores, or an equivalent AMD processor.	
RAM	16 GB	
Disk	20 GB HDD	100 GB SSD for fast loading and saving.
Graphics	No special requirements	
Operating System	Windows 10 (64 bit) Windows 11 (64 bit)	

System, Software and Configuration Compatibility

The following tables show the compatible microscope software versions, the preferred AutoScript TEM software versions per microscope software version, and the system configuration compatibility.

Although the AutoScript TEM software is backward compatible with a limited range of microscope software versions, some of the new features and improvements may be available only for the most recent supported microscope software version(s).

Compatible Microscope Software Versions per AutoScript TEM Version

It is recommended to use the highest supported version of AutoScript TEM.

AutoScript TEM	Titan	Talos	Remarks
1.0.0	3.4 - 3.7	2.4 - 2.7	
1.1.0	3.4 - 3.9	2.4 - 2.9	
1.2.0	3.6 - 3.11	2.6 - 2.11	
1.3.0	3.6 - 3.13	2.6 - 2.13	
1.4.0	3.10 - 3.15	2.10 - 2.15	
1.5.0	3.11 - 3.16	2.11 - 2.16	
1.6.0	3.12 - 3.17	2.12 - 2.17	
1.7.0	3.13 - 3.18	2.13 - 2.18	
1.8.0	3.14 - 3.19	2.14 - 2.19	
1.9.0	3.15 - 3.20	2.15 - 2.20	
1.10.0	3.16 - 3.21	2.16 - 2.21	
1.11.0	3.17 - 3.22	2.17 - 2.22	
1.12.0	3.18 - 3.23	2.18 - 2.23	
1.13.0	3.19 - 3.24	2.19 - 2.24	
1.14.0	3.20 - 3.25	2.20 - 2.25	
1.15.0	3.21 - 3.26	2.21 - 2.26	
1.16.0	3.22 - 3.27	2.22 - 2.27	

AutoScript TEM	Titan	Talos	Remarks
1.17.0	3.23 - 3.28	2.23 - 2.28	
1.18.0	3.24 - 3.29	2.24 - 2.29	

Compatible Microscope Software per AutoScript TEM Version

AutoScript TEM	Titan	Talos	Remarks
1.0.0 and later	Spectra 200	F200X	
	Spectra 300	F200S	
	Spectra Ultra	F200i	
		F200E	
		F200C	
		L120C	
1.6.0 and later	Krios G4		
	Krios G5		

Compatible Cameras for AutoScript TEM

Camera	AutoScript TEM	Remarks
Ceta 16M / -D / -M / -S	1.0.0 and later	
Ceta Speed Enhancement (Ceta-2)	1.0.0 and later	
Falcon 3EC	1.2.0 and later	
Falcon 4	1.2.0 and later	Bottom-mounted (BM) and filter-mounted (EF)

Camera	AutoScript TEM	Remarks
Falcon 4i	1.2.0 and later	Bottom-mounted (BM) and filter-mounted (EF)
EF-CCD	1.6.0 and later	

i Note: The compatible cameras may not all be supported by the compatible microscope software versions. See the TEM Server Release Notes for a definitive list of supported cameras.

Compatible Cameras and Filters for AutoScript TEM

Energy filter	AutoScript TEM	Remarks
Selectris / Selectris-X	1.3.0 and later	
Gatan Energy Filters	1.3.0 and later	Only in embedded configuration.

i Note: The compatible energy filters may not all be supported by the compatible microscope software versions. See the TEM Server Release Notes for a definitive list of supported energy filters.

Compatible Detectors for AutoScript TEM

Detector	AutoScript TEM	Remarks
HAADF	1.0.0 and later	
BF/DF Retractable	1.0.0 and later	
BF/DF Retractable Mk2	1.0.0 and later	
Panther STEM BF-S/DF-S	1.0.0 and later	
DF4	1.6.0 and later	
DualX, SuperX and UltraX	1.13.0 and later	



Note: The compatible detectors may not all be supported by the compatible microscope software versions. See the TEM Server Release Notes for a definitive list of supported detectors.

AutoScript TEM 1.18.0

New Features

Feature	Description
Access control configuration (User management)	AutoScript access control employs key-based authentication to verify the identity of clients before they are permitted to connect to the AutoScript server. Consequently, the server accepts connections exclusively from users whose public keys have been registered by an administrator. Configuring Access Control is an optional but recommended step.

Supported AutoScript TEM Objects and Functions

The following functions and objects are available in AutoScript TEM:

Static API

Functions and Objects	Supported Versions
microscope.acquisition.acquire_camera_dose_fractions()	1.16.0 and later
microscope.acquisition.acquire_camera_image()	1.0.0 and later
microscope.acquisition.acquire_camera_image_advanced()	1.2.0 and later
microscope.acquisition.acquire_camera_series()	1.0.0 and later
microscope.acquisition.acquire_camera_series_advanced()	1.2.0 and later
microscope.acquisition.acquire_stem_data	1.3.0 and later
microscope.acquisition.acquire_stem_data_advanced(settings)	1.14.0 and later
microscope.acquisition.acquire_stem_image()	1.0.0 and later
microscope.acquisition.acquire_stem_images()	1.0.0 and later
microscope.acquisition.acquire_stem_segment_images()	1.0.0 and later
microscope.acquisition.acquire_stem_segment_images_advanced()	1.8.0 and later
microscope.acquisition.continuousStemAcquisitions	1.16.0 and later
microscope.acquisition.start_acquire_stem_images_advanced()	1.16.0 and later
microscope.analysis.eds.acquire_spectrum()	1.13.0 and later
microscope.analysis.eds.acquire_spectrum_image()	1.13.0 and later
microscope.analysis.eels.acquire_spectra()	1.13.0 and later

Functions and Objects	Supported Versions
microscope.analysis.eels.acquire_spectrum_image(settings)	1.14.0 and later
microscope.auto_functions.run_auto_coma_correction()	1.2.0 and later
microscope.auto_functions.run_beam_tilt_auto_focus()	1.1.0 and later
microscope.auto_functions.run_objective_auto_stigmator()	1.1.0 and later
microscope.auto_functions.run_opti_stem()	1.4.0 and later
microscope.auto_functions.run_stem_auto_focus()	1.1.0 and later
microscope.detectors.camera_detectors	1.0.0 and later
microscope.detectors.eels_detectors	1.13.0 and later
microscope.detectors.get_camera_detector()	1.0.0 and later
microscope.detectors.get_eds_detector()	1.13.0 and later
microscope.detectors.get_eels_detector()	1.13.0 and later
microscope.detectors.get_scanning_detector()	1.0.0 and later
microscope.detectors.scanning_detectors	1.0.0 and later
microscope.detectors.screen.insert()	1.0.0 and later
microscope.detectors.screen.measure_current()	1.0.0 and later
microscope.detectors.screen.position	1.2.0 and later
microscope.detectors.screen.retract()	1.0.0 and later
microscope.optics.acceleration_voltage.limits	1.0.0 and later
microscope.optics.acceleration_voltage.value	1.0.0 and later

Functions and Objects	Supported Versions
microscope.optics.aperture_mechanisms.bi_prisms.get_available_bi_prisms()	1.12.0 and later
microscope.optics.aperture_mechanisms.bi_prisms.get_bi_prism()	1.12.0 and later
microscope.optics.aperture_mechanisms.C1	1.0.0 and later
microscope.optics.aperture_mechanisms.C2	1.0.0 and later
microscope.optics.aperture_mechanisms.C3	1.0.0 and later
microscope.optics.aperture_mechanisms.GetAvailable()	1.10.0 and later
microscope.optics.aperture_mechanisms.GetMechanism()	1.10.0 and later
microscope.optics.aperture_mechanisms.objective	1.0.0 and later
microscope.optics.aperture_mechanisms.SA	1.0.0 and later
microscope.optics.beam_diameter	1.0.0 and later
microscope.optics.beam_stopper.halfInsert()	1.12.0 and later
microscope.optics.beam_stopper.insert()	1.12.0 and later
microscope.optics.beam_stopper.insertion_state	1.12.0 and later
microscope.optics.beam_stopper.retract()	1.12.0 and later
microscope.optics.blank()	1.0.0 and later
microscope.optics.blanker.auto_beam_blank.mode.available_values	1.13.0 and later
microscope.optics.blanker.auto_beam_blank.mode.reset_to_default()	1.13.0 and later
microscope.optics.blanker.auto_beam_blank.mode.value	1.13.0 and later
microscope.optics.blanker.beam_modulator.duty_cycle	1.14.0 and later

Functions and Objects	Supported Versions
microscope.optics.blanker.beam_modulator.frequency.limits	1.14.0 and later
microscope.optics.blanker.beam_modulator.frequency.value	1.14.0 and later
microscope.optics.blanker.beam_modulator.is_modulation_on	1.14.0 and later
microscope.optics.blanker.beam_modulator.start_modulation()	1.14.0 and later
microscope.optics.blanker.beam_modulator.start_modulation(dutyCycle)	1.14.0 and later
microscope.optics.blanker.beam_modulator.start_modulation(dutyCycle, frequency)	1.14.0 and later
microscope.optics.blanker.beam_modulator.stop_modulation()	1.14.0 and later
microscope.optics.blanker.blank()	1.13.0 and later
microscope.optics.blanker.is_beam_blanked	1.13.0 and later
microscope.optics.blanker.unblank()	1.13.0 and later
microscope.optics.convergence_angle	1.0.0 and later
microscope.optics.deflectors.beam_shift	1.0.0 and later
microscope.optics.deflectors.beam_tilt	1.0.0 and later
microscope.optics.deflectors.GetAvailableDeflectors()	1.10.0 and later
microscope.optics.deflectors.GetDeflectorValue()	1.10.0 and later
microscope.optics.deflectors.gun_shift	1.0.0 and later
microscope.optics.deflectors.gun_tilt	1.0.0 and later
microscope.optics.deflectors.image_shift	1.0.0 and later

Functions and Objects	Supported Versions
microscope.optics.deflectors.image_tilt	1.0.0 and later
microscope.optics.deflectors.SetDeflectorValue()	1.10.0 and later
microscope.optics.energy_filter.apertures.available_values	1.3.0 and later
microscope.optics.energy_filter.apertures.value	1.3.0 and later
microscope.optics.energy_filter.energy_shift.limits	1.3.0 and later
microscope.optics.energy_filter.energy_shift.value	1.3.0 and later
microscope.optics.energy_filter.high_tension_offset.limits	1.3.0 and later
microscope.optics.energy_filter.high_tension_offset.value	1.3.0 and later
microscope.optics.energy_filter.mode.available_values	1.3.0 and later
microscope.optics.energy_filter.mode.prepare_mode	1.17.0 and later
microscope.optics.energy_filter.mode.value	1.3.0 and later
microscope.optics.energy_filter.slit.insert	1.3.0 and later
microscope.optics.energy_filter.slit.is_inserted	1.3.0 and later
microscope.optics.energy_filter.slit.retract	1.3.0 and later
microscope.optics.energy_filter.slit.width.limits	1.3.0 and later
microscope.optics.energy_filter.slit.width.value	1.3.0 and later
microscope.optics.focusing.stem.intensity.defocus	1.11.0 and later
microscope.optics.focusing.stem.intensity.focus	1.11.0 and later
microscope.optics.focusing.stem.intensity.reset_defocus()	1.11.0 and later

Functions and Objects	Supported Versions
microscope.optics.focusing.stem.objective.defocus	1.11.0 and later
microscope.optics.focusing.stem.objective.focus	1.11.0 and later
microscope.optics.focusing.stem.objective.reset_defocus()	1.11.0 and later
microscope.optics.focusing.tem.defocus	1.11.0 and later
microscope.optics.focusing.tem.focus	1.11.0 and later
microscope.optics.focusing.tem.reset_defocus()	1.11.0 and later
microscope.optics.illumination_mode	1.0.0 and later
microscope.optics.intensity	1.0.0 and later
microscope.optics.is_accelerator_on	1.0.0 and later
microscope.optics.is_beam_blank	1.0.0 and later
microscope.optics.is_efem_on	1.0.0 and later
microscope.optics.is_holography_mode_on	1.14.0 and later
microscope.optics.is_xl_mode_on	1.10.0 and later
microscope.optics.lenses.alignments.x	1.10.0 and later
microscope.optics.lenses.get_available()	1.9.0 and later
microscope.optics.lenses.get_lens()	1.9.0 and later
microscope.optics.magnification.available_values	1.0.0 and later
microscope.optics.magnification.increase()	1.0.0 and later
microscope.optics.magnification.value	1.0.0 and later

Functions and Objects	Supported Versions
microscope.optics.monochromator.focus	1.5.0 and later
microscope.optics.monochromator.shift	1.5.0 and later
microscope.optics.objective_lens_mode	1.0.0 and later
microscope.optics.optical_mode	1.0.0 and later
microscope.optics.paused_scan_beam_position	1.0.0 and later
microscope.optics.probe_mode	1.0.0 and later
microscope.optics.projector_mode	1.0.0 and later
microscope.optics.scan_field_of_view	1.0.0 and later
microscope.optics.scan_rotation	1.0.0 and later
microscope.optics.set_projector()	1.0.0 and later
microscope.optics.spot_size_index	1.0.0 and later
microscope.optics.stigmators.condenser_stigmator	1.1.0 and later
microscope.optics.stigmators.diffraction_stigmator	1.1.0 and later
microscope.optics.stigmators.gun_stigmator	1.1.0 and later
microscope.optics.stigmators.objective_stigmator	1.1.0 and later
microscope.optics.unblank()	1.0.0 and later
microscope.service.autoscript.client.check_pause_requested()	1.5.0 and later
microscope.service.autoscript.client.check_stop_requested()	1.5.0 and later
microscope.service.autoscript.client.is_pause_requested()	1.5.0 and later

Functions and Objects	Supported Versions
microscope.service.autoscript.client.is_stop_requested()	1.5.0 and later
microscope.service.autoscript.client.version	1.0.0 and later
microscope.service.autoscript.server.is_offline	1.3.0 and later
microscope.service.autoscript.server.version	1.0.0 and later
microscope.service.offline.reset_image	1.3.0 and later
microscope.service.offline.set_image	1.3.0 and later
microscope.service.system.name	1.0.0 and later
microscope.service.system.serial_number	1.0.0 and later
microscope.service.system.version	1.0.0 and later
microscope.source.acceleration_voltage.limits	1.16.0 and later
microscope.source.acceleration_voltage.value	1.16.0 and later
microscope.source.field_emission_gun.extractor_voltage.limits	1.16.0 and later
microscope.source.field_emission_gun.extractor_voltage.set_voltage_and_wait()	1.16.0 and later
microscope.source.field_emission_gun.extractor_voltage.value	1.16.0 and later
microscope.source.field_emission_gun.flash()	1.13.0 and later
microscope.source.field_emission_gun.gun_lens_value	1.16.0 and later
microscope.source.field_emission_gun.is_flashing_advised	1.13.0 and later
microscope.source.is_accelerator_on	1.16.0 and later
microscope.source.type	1.16.0 and later

Functions and Objects	Supported Versions
microscope.specimen.piezo_stage.disable	1.3.0 and later
microscope.specimen.piezo_stage.enable	1.3.0 and later
microscope.specimen.piezo_stage.get_axis_limits()	1.3.0 and later
microscope.specimen.piezo_stage.is_enabled	1.3.0 and later
microscope.specimen.piezo_stage.is_moving	1.3.0 and later
microscope.specimen.piezo_stage.jogging_velocity	1.3.0 and later
microscope.specimen.piezo_stage.position	1.3.0 and later
microscope.specimen.piezo_stage.relative_move()	1.3.0 and later
microscope.specimen.piezo_stage.reset_position()	1.3.0 and later
microscope.specimen.piezo_stage.start_jogging()	1.3.0 and later
microscope.specimen.piezo_stage.stop()	1.3.0 and later
microscope.specimen.stage.absolute_move()	1.3.0 and later
microscope.specimen.stage.absolute_move_safe()	1.3.0 and later
microscope.specimen.stage.get_axis_limits()	1.3.0 and later
microscope.specimen.stage.get_holder_type()	1.3.0 and later
microscope.specimen.stage.is_enabled	1.3.0 and later
microscope.specimen.stage.is_moving	1.4.0 and later
microscope.specimen.stage.move_to_pixel()	1.3.0 and later
microscope.specimen.stage.position	1.3.0 and later

Functions and Objects	Supported Versions
microscope.specimen.stage.relative_move()	1.3.0 and later
microscope.specimen.stage.relative_move_safe()	1.3.0 and later
microscope.specimen.stage.start_jogging()	1.4.0 and later
microscope.specimen.stage.stop_jogging()	1.4.0 and later
microscope.stage.absolute_move()	1.0.0 - 1.3.0
microscope.stage.absolute_move_safe()	1.0.0 - 1.3.0
microscope.stage.get_axis_limits()	1.0.0 - 1.3.0
microscope.stage.get_holder_type()	1.0.0 - 1.3.0
microscope.stage.is_enabled	1.0.0 - 1.3.0
microscope.stage.position	1.0.0 - 1.3.0
microscope.stage.relative_move()	1.0.0 - 1.3.0
microscope.stage.relative_move_safe()	1.0.0 - 1.3.0
microscope.vacuum.column_valves.close()	1.2.0 and later
microscope.vacuum.column_valves.open()	1.2.0 and later
microscope.vacuum.column_valves.state	1.2.0 and later
microscope.vacuum.state	1.0.0 and later

Dynamic Objects

Functions and Objects	Supported Versions
ApertureMechanism.aperture	1.0.0 and later

Functions and Objects	Supported Versions
ApertureMechanism.apertures	1.0.0 and later
ApertureMechanism.disable()	1.0.0 and later
ApertureMechanism.enable()	1.0.0 and later
ApertureMechanism.insert()	1.0.0 and later
ApertureMechanism.insertion_state	1.4.0 and later
ApertureMechanism.is_enabled	1.0.0 and later
ApertureMechanism.is_retractable	1.0.0 and later
ApertureMechanism.position	1.0.0 and later
ApertureMechanism.reset_positions()	1.0.0 and later
ApertureMechanism.retract()	1.0.0 and later
BiPrism.disable()	1.12.0 - 1.13.0
BiPrism.enable()	1.12.0 - 1.13.0
BiPrism.is_enabled	1.12.0 and later
BiPrism.name	1.12.0 and later
BiPrism.recover()	1.12.0 and later
BiPrism.rotation	1.12.0 and later
BiPrism.voltage	1.12.0 and later
BiPrismRotation.is_available	1.12.0 and later
BiPrismRotation.limits	1.12.0 and later
BiPrismRotation.rotate(angle)	1.12.0 and later

Functions and Objects	Supported Versions
BiPrismRotation.Value	1.12.0 and later
BiPrismVoltage.limits	1.12.0 and later
BiPrismVoltage.value	1.12.0 and later
CameraDetector.display_name	1.1.0 and later
CameraDetector.exposure_time_limits	1.1.0 and later
CameraDetector.get_supported_image_sizes()	1.14.0 and later
CameraDetector.insert()	1.1.0 and later
CameraDetector.insertion_state	1.1.0 and later
CameraDetector.is_dose_sensitive	1.1.0 and later
CameraDetector.is_insertable	1.1.0 and later
CameraDetector.is_operational	1.1.0 and later
CameraDetector.name	1.1.0 and later
CameraDetector.pixel_size	1.1.0 and later
CameraDetector.retract()	1.1.0 and later
ContinuousStemAcquisition.ContinueAcquisition()	1.16.0 and later
ContinuousStemAcquisition.GetAllImages()	1.16.0 and later
ContinuousStemAcquisition.GetImage()	1.16.0 and later
ContinuousStemAcquisition.State	1.16.0 and later
ContinuousStemAcquisition.StopAcquisition()	1.16.0 and later
ContinuousStemAcquisition.WaitForNextFrame()	1.16.0 and later

Functions and Objects	Supported Versions
DetectorSegment.detector_name	1.0.0 and later
DetectorSegment.gain	1.0.0 and later
DetectorSegment.is_enabled	1.0.0 and later
DetectorSegment.name	1.0.0 and later
DetectorSegment.offset	1.0.0 and later
Lens.limits_raw	1.9.0 and later
Lens.name	1.9.0 and later
Lens.value	1.9.0 and later
Lens.value_raw	1.9.0 and later
PresetRaw.limits	1.10.0 and later
PresetRaw.value	1.10.0 and later
ScanningDetector.collection_angle_limits	1.0.0 and later
ScanningDetector.display_name	1.1.0 and later
ScanningDetector.gain	1.0.0 and later
ScanningDetector.gain_limits	1.6.0 and later
ScanningDetector.get_enabled_segments()	1.0.0 and later
ScanningDetector.get_segment()	1.0.0 and later
ScanningDetector.get_segments()	1.0.0 and later
ScanningDetector.insert()	1.1.0 and later
ScanningDetector.insertion_state	1.1.0 and later

Functions and Objects	Supported Versions
ScanningDetector.is_insertable	1.6.0 and later
ScanningDetector.name	1.0.0 and later
ScanningDetector.offset	1.0.0 and later
ScanningDetector.offset_limits	1.6.0 and later
ScanningDetector.relative_gain	1.6.0 and later
ScanningDetector.relative_offset	1.6.0 and later
ScanningDetector.reset_segments()	1.0.0 and later
ScanningDetector.retract()	1.1.0 and later
ScanningDetector.set_enabled_segments()	1.0.0 and later
ScanningDetector.set_segment_gain()	1.0.0 and later
ScanningDetector.set_segment_offset()	1.0.0 and later
XLensAlignment.preset_raw	1.10.0 and later

Other Product Information

AutoScript 1 is available in 3 different product variants for the Python scripting usecase. Additionally, there is one feature variant that can be combined with the product variant. The variants determine which microscope parameters are accessible for scripting. All variants are installed with the same installation media, the distinction is made by applying a specific license key after installation. The variants are:

Product Variant Name	Commercial ID
AutoScript TEM 1 - Academic	1408417
AutoScript TEM 1 - Commercial	1408419
AutoScript TEM 1 - Offline	1290580

Feature Variant Name	Commercial ID

It is possible to evaluate the product during the 90-day evaluation period without a licence key. For the duration of this trial period the "AutoScript TEM 1 - Commercial " variant is activated.

Only one product variant can be activated on a single computer with combination of any number of feature product variants.

Known Issues

The most recent release of AutoScript TEM has the following Known Issues.

ID	Issue Description and Workaround
ASTEM-2169	New environment cannot be created.
	Run the PackageManager as administrator.

ID	Issue Description and Workaround
ASTEM-1972	SmartCam viewer in Velox is blocking Flucam acquisitions.
	Stop the SmartCam preview manually.

ID	Issue Description and Workaround
ASTEM-1861	Bi-prism interface is not synchronous.
	Add delays between calls.

ID	Issue Description and Workaround
ASTEM-1543	Move to pixel does not work properly with non-zero STEM rotation.
	Do not use STEM scan rotation.